Español

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English

Introduction

Congratulations for having purchased a Drive Medical Ultrasonic Nebulizer!  
This ultrasonic nebulizer is intended to spray liquids in aerosol form into gases that are delivered directly to the patient for breathing for use by both the adult and pediatric populations.  
Drive Medical warrants that this device has been produced using high quality materials and it complies with the safety rules in force both domestically and internationally.  
The Drive Medical ultrasonic nebulizer **operates on ultrasound which produces the drug solution vibration and turns it into mist, these tiny particles get deposited in the lungs and consequently, effective and quick results are achieved and no side- effects are caused.**  
This model comes with the **exclusive Closed Water Chamber System**, this system avoids filling with water each time. It is easier to use: **: all you need to do is to pour the solution into the medication cup and it is ready to use.**  
Silent and safe, it stops working automatically 3 minutes after the unit was switched on, preventing the machine from working without medication.  
It weighs only 57 lbs (approx.), it has been designed to operate with different power supplies which allows use in different situations: external source, apt to be connected to power supply at home, rechargeable battery pack with rechargeable batteries and car lighter adapter (optional).

**CAUTION: As per USA Federal regulations this product can only be sold under medical prescription.**

Important Safeguards

- 1- Read the manual carefully before using the device and keep it in case you need further information in the future.**
- 2- This product was intended for aerosoltherapy only. Any other use will be considered wrong and consequently dangerous; Drive Medical will not be held responsible for damages caused by misuse or if the device has been used in electric fittings which do not comply with the safety rules in force.**
- 3- Do not turn the device on if there is no drug solution into the medication cup.**
- 4- Disconnect the nebulizer when it is not being used.**
- 5- Do not remove the medication cup cover while the device is operating.**
- 6- Do not expose the nebulizer to a heat source.**
- 7- Do not pour more than 8 ml. of solution into the medication cup.**
- 8- Do not leave drug solution in the medication cup when the device is not being used.**
- 9- Do not scrape or pierce the medication cup cover.**
- 10- Do not touch the medication cup or the drug solution while the device is being used.**
- 11- Do not dip the nebulizer in water.**
- 12- Children may only use the device under adult supervision.**
- 13- To avoid infections and contaminations the accessories should not be shared by different patients. They are intended for personal use.**
- 14- Use the unit intermittently. Maximum operating period: 30 minutes. Rest period: 30 minutes.**
- 15- This device can only be repaired by an authorized Drive dealer.**
- 16- Never operate this nebulizer if it has a damaged cord or plug, if it is NOT working properly, if it has been dropped or damaged, or dropped into water. Call qualified Service Personnel for examination and repair.**

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Accessories Included

1 Pediatric mask / 1 Adult mask / 1 Adapter / 1 Nasal piece / 2 Air Filters / 1 Mouth piece / 1 AC Adapter / Rechargeable Battery Pack / 1 Instructions Manual and Warranty

Optional Accessories: car lighter adapter

Spare parts descriptions

See pictures on page 2.

Use Instructions

(See pictures on page 3)

- 1- Take the unit and the AC adapter out of the box
  - 2- Remove the clear housing and pour the medication into the medication container (8 ml. Max.)
  - 3- Put the clear housing and insert the device to be used into it (adult mask, pediatric mask, and nasal or mouth piece).
- CAUTION:** carefully check the membrane placed at the bottom of the medication cup before each use. If it is pierced or broken do not use the nebulizer and contact the authorized Service Agent.
- 4- Plug the unit and the power connector to the unit.
  - 5- Press the ON switch, the green light will turn on and the unit will start to operate. Adjust the air regulator to the desired intensity.
  - 6- The unit will automatically turn off at 3 minutes. If you wish to continue nebulizing, check if the medication amount in the pool is the suitable one. Add more if necessary. Then press the ON/OFF switch twice and the unit will start to operate again.
  - 7- Once use is over, unplug the equipment and disconnect the AC adapter. Remove the clear housing and dry the left medication using a clean cloth.

Operation with Car Lighter Adapter (Optional accessory)

Connect the adapter to the car lighter and the other wire end to the unit. Then operate the nebulizer as per the instructions previously described.

Operation with Battery Pack

Suggestions and Safety measures

- 1- Do not use in direct sunlight or near heat source.
- 2- Never submerge the Battery Pack in water.
- 3- Do not use the battery pack for more than 20 consecutive minutes.
- 4- To start, be sure battery pack is fully charged, (approx. 14 hrs).
- 5- Plug the battery pack into the ultrasonic port. Switch the power pack on "I".
- 6- When complete, turn off battery pack, remove plug.
- 7- If you did not use the battery pack for 3 months or more, recharge.
- 8- Caution: the battery pack temperature may increase while in use.

English

English

Initial charge of the Rechargeable Battery Pack

The Pack is provided uncharged and it is necessary to be completely charge it before using it for the first time to ensure the capacity and the life spun of the batteries.

- 1- Make sure that the start button is in OFF “0” position.
- 2- Connect the output of your compact nebulizer to the mains.
- 3- Connect the plug of the output in the hole placed on the bottom of the pack and leave it charging for 14 hours.

Regular use charge

- 1- Make sure that the start button is in OFF position “0”.
- 2- Connect the output of your compact nebulizer to the mains.
- 3- Connect the plug of the output in the hole placed on the bottom of the pack and leave it charging between 6 and 12 hours, depending on the use. Do not leave the charge connected for more than 24 hours.

Instructions for use

Prepare the nebulizer as you usually do or follow the corresponding instructions. Then connect the plug of the pack to the nebulizer and move the on-off button to ON position “I”, the green led will turn on to indicate the condition of the batteries. Turn on the nebulizer by pushing the start button. If the red light turns on, disconnect the nebulizer plug and connect the pack to be charged. With the complete charge, the pack autonomy is of approx. 20 minutes. When you finish using the nebulizer, turn off the pack and disconnect the plug.

Technical Specifications

Type of batteries .....	Nickel – metal hydride
Capacity .....	550mAh
Output .....	24V $\Rightarrow$ 0,5A
Input .....	20V $\Rightarrow$ 0,5A
Charge time .....	14 hs first use 6 to 12 hs regulare use
Operating conditions .....	(50°F~104°F) +10°C~+40°C/ 45~85 RH%
Storage conditions .....	(-4°F~113°F) -20°C~+45°C/ 45~85 RH%
Weight .....	8 oz

The batteries material must be recycled or deposited in a proper place to such end.

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Clean and Maintenance

**Important:** Clean the accessories before using the nebulizer for the first time or after each use following these instructions: wash all the pieces, except for the air filter, with warm water and detergent. Rinse with warm tap water to remove the detergent. Dip the pieces into a clean bowl containing 3 parts of warm water and 1 part of vinegar during 30 minutes, disinfectant can also be used. Wash your hands and take out the pieces from the disinfectant solution, rinse them with warm tap water and leave them on an absorbent piece of paper. Keep the pieces into a hermetic bag. Do not dry them with a towel. Clean the outer surface with a dry piece of cloth. CAUTION: Do not boil. **WARNING: In order to avoid infections, it is highly recommended to clean the pieces after each use. Only nebulize solutions prescribed by a doctor and follow his instructions carefully. Do not use balm substances. Do not use any volatile fluid as alcohol or benzene** The life time of the accessories, using them according to the instructions of this manual is of 5 years. 1- The air filter has to be replaced if it turns dirty or torn. 2- Remove the air chamber from the medication cup cover. 3- Take out the filter and place a new one. 4- Put all the pieces back .

Technical Specifications

Method of nebulization.....	Ultrasonic Oscillation
Operating voltage.....	12V
Class.....	
Power.....	18 W
Oscillation frequency.....	2.5 Mhz
Nebulization rate.....	0.5 ml/ min(minimum)**
Particle dimension.....	1.5 to 5.7 µm (MMAD 3.8 µm)
Medication cup maximum capacity.....	8 ml
Weight.....	0.57 lb / 260 g
Device dimensions.....	3.2" x 6.6" x 2.3" 83.5 x 169 x 60 mm
Sound noise level (at 1m distance).....	.25 dB
Operating temperature and humidity.....	41 to 104°F, 30% to 95% 5 to 40° C, 30% to 95%
Storage temperature and humidity.....	33.8 to 104°F, 10% to 95% 1 to 40° C, 10% to 95%
Atmospheric pressure range for transport and storage.....	50 to 106 kPa (500-1060 mbar)
Service Cycle.....	30 minutes on, 30 minutes off
Noise level (at 1 meter).....	.25 dB

\*The Closed Water Chamber contains a solution of Benzalkonium Chloride at 0,01%. \*The specifications and the outer design are subject to change for improvement reasons without previous notice. \*\* Value measured with 0,9% of ClNa solution. This specification may change based on the type and viscosity of the medication used.

**The Compact Ultrasonic Nebulizer Model 3060 operates with the external power supply provided by San Up:** Primary voltage..... 100-240V 50/60Hz 0.45 A Secondary voltage..... 12V  $\approx$  o 1,25A

English

### Spare parts list

See on page 4. Please indicate the parts code number when ordering spare parts.

### Trouble Shooting Chart

In the event the unit starts malfunctioning or stops, review the following:

Symptom	Possible cause	Suggested solutions
Although the device is on, it does not work	<b>1.</b> Not plugged in. <b>2.</b> No power <b>3.</b> The AC Adapter is not connected	<b>1.</b> Check the plug is properly connected <b>2.</b> Check the power is not gone out <b>3.</b> Check the AC Adapter is properly connected
The mist does not come out	There is no solution into the cup	Pour some drug solution into the medication cup
Little mist is delivered	<b>1.</b> The air filter is dirty <b>2.</b> There is no enough solution into the cup <b>3.</b> Little power <b>4.</b> Air excess	<b>1.</b> Change the filter <b>2.</b> Pour some drug solution into the medication cup <b>3.</b> Increase the power <b>4.</b> Regulate the air flow
The device operates for a period and then turns off	Timer off	Press the key twice and turn it on again

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Certifications

This product meets the highest international quality standards.

ISO 9001 Quality Systems (International)

13485 Medical Equipment (European Community)

A.N.M.A.T. (Argentina)

(Argentina)

Manufactured for Drive Medical.

Made in Argentina

English

English

### 3 Year Limited Warranty

Drive Medical warrants to the original retail purchaser that this product is free from defects in material and workmanship for a period of 3 years from the date of purchase. This warranty does not cover failures resulting from misuse, abuse, alteration of the Drive Medical Nebulizer System, use improper components or failure to comply with the operating instructions. This warranty does not include the Nebulizer System tubing or accessories.

THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER ORAL, WRITTEN, EXPRESS OR IMPLIED. THERE ARE NO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, EXCEPT FOR THOSE IMPLIED WARRANTIES ARISING BY OPERATION OF LAW WHICH ARE LIMITED IN DURATION TO THE 3 YEAR TERM OF THIS WARRANTY. DRIVE MEDICAL'S LIABILITY, WHETHER BASED ON CONTRACT, TORT, WARRANTY, STRICT LIABILITY, OR ANY OTHER THEORY, SHALL NOT EXCEED THE PRICE OF THE INDIVIDUAL DRIVE MEDICAL NEBULIZER SYSTEM WHICH IS THE BASIS OF THE CLAIM. IN NO EVENT SHALL DRIVE MEDICAL BE LIABLE FOR ANY LOSS OF PROFITS, FACILITIES, EQUIPMENT OR OTHER INCIDENTAL, CONSEQUENTIAL OR SPECIAL DAMAGES.

Some states do not allow limitations on how long an implied warranty lasts, or allow the exclusion or limitation of incidental, consequential, or special damages, so the above limitations or exclusions may not apply to you. This warranty gives you specific legal rights. You may also have other rights which vary from state to state.

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